

Dengue 2 NS1 Antibody

Polyclonal Rabbit Anti Dengue 2 NS1
ABP0007

Product Overview

Name Dengue 2 NS1 Antibody

Description

Polyclonal Rabbit Anti Dengue 2 NS1

Source

Rabbit.

Formulation

200mM glycine, pH 7.5 and 0.02% NaN₃.

Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Purity

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Purification Method

Purified by protein A chromatography.

Precautions

Dengue 2 NS1 Antibody is for research use only and not for use in diagnostic or therapeutic procedures.