

DiMethyl-Histone H3 (Lys36) Rabbit mAb


ABZ4618

Product Overview

Name	DiMethyl-Histone H3 (Lys36) Rabbit mAb
Catalog #	ABZ4618
Clonality	Monoclonal
Accession(Primary)	P6843
Application Note (Approx.)	WB1:1000-1:5000
Precautions	

Target information(P6843)

Synonyms**Gene ID** 8350**Other Names****Function****Cellular location****Note**



Western blot analysis of DiMethyl-Histone H3 (Lys36) in HeLa, 3T3 lysates using DiMethyl-Histone H3 (Lys36) antibody.