

# DiMethyl-Histone H3 (Lys4) Rabbit mAb


ABZ5297

## Product Overview


<b>Name</b>	DiMethyl-Histone H3 (Lys4) Rabbit mAb
<b>Catalog #</b>	ABZ5297
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P6843
<b>Application Note (Approx.)</b>	WB1:500-1:1000 IHCP-1:200-1:1000
<b>Precautions</b>	

## Target information(P6843)

**Synonyms****Gene ID** 8350**Other Names****Function****Cellular location****Note**



Western blot analysis of MethylHistone H3 (di K4) in HeLa lysates using DiMethyl-Histone H3 (Lys4) antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using MethylHistone H3 (di K4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.