

# DiMethyl-Histone H3 (Lys79) (8E8) Mouse mAb

ABZ4603

## Product Overview

<b>Name</b>	DiMethyl-Histone H3 (Lys79) (8E8) Mouse mAb
<b>Catalog #</b>	ABZ4603
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P6843
<b>Application Note (Approx.)</b>	WB1:500-1:1000
<b>Precautions</b>	

## Target information(P6843)

**Synonyms****Gene ID** 8350**Other Names****Function****Cellular location****Note**



Western blot analysis of DiMethyl-Histone H3 (Lys79) in HeLa lysates using DiMethyl-Histone H3 (Lys79) antibody.