

DiMethyl-Histone H3 (Lys9) Rabbit mAb

ABZ4622

Product Overview

Name	DiMethyl-Histone H3 (Lys9) Rabbit mAb
Catalog #	ABZ4622
Clonality	Monoclonal
Accession(Primary)	P68431
Application Note (Approx.)	WB1:500-1:1000
Precautions	

Target information(P68431)

Synonyms**Gene ID** 8350**Other Names****Function****Cellular location****Note**



Western blot analysis of DiMethyl-Histone H3 (Lys9) in HeLa, lysates using DiMethyl-Histone H3 (Lys9) antibody.