

## H5N1 Polyclonal Antibody

*Influenza-A Hemagglutinin H5N1, Polyclonal Rabbit Antibody  
ABP0011*

### Product Overview

Name H5N1 Polyclonal Antibody

#### Description

Influenza-A Hemagglutinin H5N1, Polyclonal Rabbit Antibody

#### Source

Baculovirus Insect Cells.

#### Physical Appearance

Sterile Filtered colorless solution.

#### Formulation

The Recombinant H7N7 A/Netherlands/219/03 solution contains 10mM Sodium phosphate, pH 7 and 150mM NaCl.

#### Stability

H7N7 A/Netherlands/219/03 Recombinant should be stored at 4°C.

#### Purity

Greater than 90.0% as determined by SDS-PAGE.

#### Immunological Activity

Western-Blot 0.1 µg -1 µg per strip, ELISA 1 µg/Well.

#### Precautions

H5N1 Polyclonal Antibody is for research use only and not for use in diagnostic or therapeutic procedures.