

# HA Antibody

HA Polyclonal Antibody  
ABP0012

## Product Overview

Name HA Antibody

### Description

HA Polyclonal Antibody

### Physical Appearance

Sterile filtered colorless solution.

### Formulation

1xPBS & 50% glycerol.

### Applications

Western Blot, Immunoprecipitation.

### Type

Mouse Antibody Monoclonal.

### Clone

PHASHG.

### Ig Subclass

Mouse IgG1.

### Storage Procedures

For periods up to 1 month store at 4°C, for longer periods of time, store at -20°C. Prevent freeze thaw cycles.

### Titer

Western Blotting: 1 µg/ml.

### Precautions

HA Antibody is for research use only and not for use in diagnostic or therapeutic procedures.