

# KD-Validated Anti-ABCD1 Rabbit Monoclonal Antibody

*Rabbit monoclonal antibody*

ABG1172

## Product Overview

Name	KD-Validated Anti-ABCD1 Rabbit Monoclonal Antibody
Catalog #	ABG1172
Clonality	Monoclonal
Accession(Primary)	P33897
Application Note (Approx.)	WB1:5,000
Precautions	

## Target information(P33897)

Synonyms


Gene ID

Other Names


Function

Cellular location

Note



Western blotting analysis using anti-ABCD1 antibody (Cat#ABG1172). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-ABCD1 antibody (Cat#ABG1172, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-ABCD1 antibody (Cat#ABG1172). ABCD1 expression in wild type (WT) and ABCD1 shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-ABCD1 antibody (Cat#ABG1172, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.