

KD-Validated Anti-ABCF2 Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG1955

Product Overview

Name	KD-Validated Anti-ABCF2 Mouse Monoclonal Antibody
Catalog #	ABG1955
Clonality	Monoclonal
Accession(Primary)	Q9UG63
Application Note (Approx.)	WB1:5,000
Precautions	

Target information(Q9UG63)

Synonyms


Gene ID

Other Names

Function

Cellular location

Note





Western blotting analysis using anti-ABCF2 antibody (Cat#ABG1955). ABCF2 expression in wild type (WT) and ABCF2 shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-ABCF2 antibody (Cat#ABG1955, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.