

KD-Validated Anti-ACACA Rabbit Polyclonal Antibody

Rabbit polyclonal antibody

ABG2142

Product Overview

Name	KD-Validated Anti-ACACA Rabbit Polyclonal Antibody
Catalog #	ABG2142
Clonality	Polyclonal
Accession(Primary)	Q13085
Application Note (Approx.)	WB1:2,500
Precautions	

Target information(Q13085)

Synonyms


Gene ID

Other Names


Function

Cellular location

Note



Western blotting analysis using anti-ACACA antibody (Cat#ABG2142). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-ACACA antibody (Cat#ABG2142, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-ACACA antibody (Cat#ABG2142). ACACA expression in wild-type (WT) and ACACA shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 serves as a loading control. The blot was incubated with anti-ACACA antibody (Cat#ABG2142, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody respectively.