

KD-Validated Anti-ACAT1 Rabbit Polyclonal Antibody

Rabbit polyclonal antibody


ABG2113

Product Overview


Name	KD-Validated Anti-ACAT1 Rabbit Polyclonal Antibody
Catalog #	ABG2113
Clonality	Polyclonal
Accession(Primary)	P24752
Application Note (Approx.)	WB1:2,500
Precautions	

Target information(P24752)

Synonyms**Gene ID****Other Names****Function****Cellular location****Note**



Western blotting analysis using anti-ACAT1 antibody (Cat#ABG2113). Total cell lysates (20 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-ACAT1 antibody (Cat#ABG2113, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-ACAT1 antibody (Cat #ABG2113). ACAT1 expression in wild-type (WT) and ACAT1 shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. Hsp90 serves as a loading control. The blot was incubated with anti-ACAT1 antibody (Cat #ABG2113, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody respectively.