

KD-Validated Anti-Basigin Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1799

Product Overview


Name	KD-Validated Anti-Basigin Rabbit Monoclonal Antibody
Catalog #	ABG1799
Clonality	Monoclonal
Accession(Primary)	P35613
Application Note (Approx.)	WB1:5,000 FC1:2,000 ICC1:1,000
Precautions	

Target information(P35613)


Synonyms**Gene ID****Other Names****Function****Cellular location****Note**




Western blotting analysis using anti-CD147 antibody (Cat#ABG1799). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-CD147 antibody (Cat#ABG1799, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-CD147 antibody (Cat#ABG1799). CD147 expression in wild-type (WT) and BSG shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 serves as a loading control. The blot was incubated with anti-CD147 antibody (Cat#ABG1799, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of CD147 expression in HepG2 cells using anti-CD147 antibody (Cat#ABG1799, 1:2,000). Green, isotype control; red, CD147.



Immunocytochemical staining of HepG2 cells with anti-CD147 antibody(Cat#ABG1799, 1:1,000). Nuclei were stained blue with DAPI; CD147 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.
