

KD-Validated Anti-active + pro caspase 3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2292

Product Overview


Name	KD-Validated Anti-active + pro caspase 3 Rabbit Monoclonal Antibody
Catalog #	ABG2292
Clonality	Monoclonal
Accession(Primary)	P42574
Application Note (Approx.)	WB1:5,000
Precautions	

Target information(P42574)

Synonyms**Gene ID****Other Names****Function****Cellular location****Note**



Western blotting analysis using anti-pro Caspase 3 antibody (Cat#ABG2292). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-pro Caspase 3 antibody (Cat#ABG2292, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-pro Caspase 3 antibody (Cat#ABG2292). pro Caspase 3 expression in wild type (WT) and pro Caspase 3 shRNA knockdown (KD) HeLa cells with 30 μ g of total cell lysates. β -Tubulin serves as a loading control. The blot was incubated with anti-pro Caspase 3 antibody (Cat#ABG2292, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.