

# alpha Smooth Muscle Actin Rabbit pAb


ABZ1369

## Product Overview

<b>Name</b>	alpha Smooth Muscle Actin Rabbit pAb
<b>Catalog #</b>	ABZ1369
<b>Clonality</b>	Polyclonal
<b>Accession(Primary)</b>	P62736
<b>Application Note (Approx.)</b>	WB1:500-1:2000 IHC1:100-1:300
<b>Precautions</b>	

## Target information(P62736)


**Synonyms****Gene ID** 70**Other Names****Function****Cellular location****Note**



Western blot analysis of ACTC1 in MFC lysates using ACTC1 antibody.



Immunohistochemistry analysis of paraffin-embedded rat muscle using ACTC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded rat muscle using ACTC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse muscle using ACTC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse muscle using ACTC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.